

TMPRSS2 antibody



Catalog Number: 116274

Product name

TMPRSS2 antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat; other species not tested.

Applications

ELISA, WB

Antibody description

TMPRSS2 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells. Observed molecular weight by Western-blot: 54 kDa

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of TMPRSS2 recombinant protein (Accession Number: XM_011529731). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:116274(TMPRSS2 antibody) at dilution of 1:1000

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality