

TMTC1 antibody



Catalog Number: 116281

Product name

TMTC1 antibody

Immunogen

[Human TMTC1 Recombinant protein \(His tag\)](#)

Specificity

human, mouse; other species not tested.

Antibody description

TMTC1 Rabbit Polyclonal antibody. Positive WB detected in mouse bladder tissue, mouse spleen tissue, mouse colon tissue. Observed molecular weight by Western-blot: 99 kDa, 65 kDa

Preparation

This antibody was obtained by immunization of TMTC1 recombinant protein (Accession Number: NM_175861). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



mouse bladder tissue were subjected to SDS PAGE followed by western blot with Catalog No:116281(TMTC1 antibody) at dilution of 1:1000