

TNFAIP3 antibody

Catalog Number: 116286

Product name

TNFAIP3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TNFAIP3 Rabbit Polyclonal antibody. Positive WB detected in RAW 264.7 cells, HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 90 kDa, 37-40 kDa

Preparation

This antibody was obtained by immunization of TNFAIP3 recombinant protein (Accession Number: NM_001270507). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



RAW264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:116286(TNFAIP3 antibody) at dilution of 1:1000