

TNFR1 antibody

Catalog Number: 116289

Product name

TNFR1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TNFR1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HL-60 cells, human brain tissue, NIH/3T3 cells, Raji cells. Positive FC detected in Raji cells. Observed molecular weight by Western-blot: 40-45 kDa, 60 kDa

Preparation

This antibody was obtained by immunization of TNFR1 recombinant protein (Accession Number: NM_001065). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, FC

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116289(TNFR1 antibody) at dilution of 1:500



1X10⁶ Raji cells were stained with 0.2ug TNFR1 antibody (Catalog No:116289, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.