

TNFSF13B antibody



Catalog Number: 116294

Product name

TNFSF13B antibody

Specificity

human, mouse; other species not tested.

Antibody description

TNFSF13B Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, human liver tissue, U-937 cells. Observed molecular weight by Western-blot: 32kd,45-52kd

Preparation

This antibody was obtained by immunization of TNFSF13B recombinant protein (Accession Number: NM_006573). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116294(TNFSF13B antibody) at dilution of 1:1000