

TRIM33 antibody



Catalog Number: 116312

Product name

TRIM33 antibody

Specificity

Human; other species not tested.

Antibody description

TRIM33 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells, MCF7 cells, PC-3 cells. Observed molecular weight by Western-blot: 140-150 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_015906). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:116312(TRIM33 antibody) at dilution of 1:500