

# TRIM38 antibody



Catalog Number: 116314

## Product name

TRIM38 antibody

## Specificity

Human; other species not tested.

## Antibody description

TRIM38 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, HEK-293 cells, human colon tissue. Observed molecular weight by Western-blot: 53 kDa

## Preparation

This antibody was obtained by immunization of TRIM38 recombinant protein (Accession Number: XM\_024446303). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:116314(TRIM38 antibody) at dilution of 1:800