

# TRIM5 antibody



Catalog Number: 116319

## Product name

TRIM5 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

TRIM5 Rabbit Polyclonal antibody. Positive WB detected in MKN-45 cells, mouse stomach tissue. Positive IHC detected in human prostate cancer tissue, human stomach tissue. Observed molecular weight by Western-blot: 46 kDa, 70 kDa

## Preparation

This antibody was obtained by immunization of TRIM5 recombinant protein (Accession Number: NM\_033092). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:116319(TRIM5 antibody) at dilution of 1:50



(under 10x lens)

MKN-45 cells were subjected to SDS PAGE followed by western blot with Catalog No:116319(TRIM5 Antibody) at dilution of 1:300