

TRIM69 antibody



Catalog Number: 116324

Product name

TRIM69 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TRIM69 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue. Positive WB detected in mouse testis tissue, HeLa cells. Observed molecular weight by Western-blot: 30kd

Preparation

This antibody was obtained by immunization of TRIM69 recombinant protein (Accession Number: NM_080745). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116324(TRIM69 Antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:116324(TRIM69 antibody) at dilution of 1:100 (under 10x lens)