

TRAPPC2 antibody

Catalog Number: 116342

Product name

TRAPPC2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TRAPPC2 Rabbit Polyclonal antibody. Positive IP detected in PC-3 cells. Positive WB detected in human brain tissue, human liver tissue, human lung tissue, mouse brain tissue, mouse thymus tissue, PC-3 cells. Observed molecular weight by Western-blot: 16 kDa

Preparation

This antibody was obtained by immunization of TRAPPC2 recombinant protein (Accession Number: XM_011545565). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116342 (TRAPPC2 antibody) at dilution of



1:600

IP Result of anti-TRAPPC2 (IP:Catalog No:116342, 3ug; Detection:Catalog No:116342 1:500) with PC-3 cells lysate 2800ug.