

TRAPPC3 antibody



Catalog Number: 116344

Product name

TRAPPC3 antibody

Immunogen

[Human TRAPPC3 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

TRAPPC3 Rabbit Polyclonal antibody. Positive IP detected in mouse liver tissue. Positive WB detected in mouse liver tissue, HEK-293 cells, HepG2 cells, mouse small intestine tissue, PC-3 cells. Positive IHC detected in human placenta tissue. Positive IF detected in Hela cells. Observed molecular weight by Western-blot: 20-22 kDa

Preparation

This antibody was obtained by immunization of TRAPPC3 recombinant protein (Accession Number: NM_014408). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF, IHC, IP

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000


IHC: 1:20-1:200

IF: 1:10-1:100

Validations




mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:116344(TRAPPC3 antibody) at dilution of




1:800
Immunohistochemical of paraffin-embedded human placenta using Catalog No:116344(TRAPPC3 antibody) at dilution of 1:100 (under 10x lens)

TRAPPC3 antibody

Catalog Number: 116344



Immunofluorescent analysis of HeLa cells, using TRAPPC3 antibody Catalog No:116344 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

IP Result of anti-TRAPPC3 (IP:Catalog No:116344, 3ug; Detection:Catalog No:116344 1:800) with mouse liver tissue lysate 4000ug.