

TSPAN31 antibody



Catalog Number: 116361

Product name

TSPAN31 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TSPAN31 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells, mouse skeletal muscle tissue, PC-3 cells. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of TSPAN31 recombinant protein (Accession Number: NM_005981). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116361(TSPAN31 antibody) at dilution of 1:400