

TSPY3 antibody

Catalog Number: 116367

Product name

TSPY3 antibody

Immunogen

[Human TSPY3 Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

TSPY3 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue. Positive WB detected in human testis tissue, LNCaP cells. Observed molecular weight by Western-blot: 39-40 kDa

Preparation

This antibody was obtained by immunization of TSPY3 recombinant protein (Accession Number: NM_001077697). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116367(TSPY3 antibody) at dilution of 1:2000



Immunohistochemical of paraffin-embedded human testis using Catalog No:116367(TSPY3 antibody) at dilution of 1:100 (under 10x lens)