

TTC26 antibody



Catalog Number: 116377

Product name

TTC26 antibody

Immunogen

[Human TTC26 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

TTC26 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, HEK-293 cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 60-64 kDa

Preparation

This antibody was obtained by immunization of TTC26 recombinant protein (Accession Number: NM_024926). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116377(TTC26 Antibody) at dilution of 1:300