

TTC30B antibody



Catalog Number: 116379

Product name

TTC30B antibody

Immunogen

[Human TTC30B Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

TTC30B Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, mouse heart tissue. Observed molecular weight by Western-blot: 68-70 kDa

Preparation

This antibody was obtained by immunization of TTC30B recombinant protein (Accession Number: NM_152517). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116379(TTC30B Antibody) at dilution of 1:600