

TTC32 antibody

Catalog Number: 116381

Product name

TTC32 antibody

Immunogen

[Human TTC32 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

TTC32 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Positive IP detected in mouse testis tissue. Positive IF detected in HEK-293 cells. Positive IHC detected in human testis tissue. Observed molecular weight by Western-blot: 17-30 kDa

Preparation

This antibody was obtained by immunization of TTC32 recombinant protein (Accession Number: NM_001008237). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP, IF

Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000


IHC: 1:20-1:200

IF: 1:10-1:100

Validations




mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116381(TTC32 Antibody) at dilution of 1:300




IP Result of anti-TTC32 (IP:Catalog No:116381, 4ug; Detection:Catalog No:116381 1:300) with mouse testis tissue lysate 4000ug.

TTC32 antibody

Catalog Number: 116381



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:116381(TTC32 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HEK-293 cells using Catalog No:116381(TTC32 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)