

TTC35 antibody

Catalog Number: 116382

Product name

TTC35 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TTC35 Rabbit Polyclonal antibody. Positive IP detected in mouse heart tissue. Positive WB detected in Jurkat cells, human placenta tissue, K-562 cells, MCF-7 cells, mouse heart tissue, mouse spleen tissue, NIH/3T3 cells, PC-3 cells, SKOV-3 cells. Observed molecular weight by Western-blot: 35 kDa

Preparation

This antibody was obtained by immunization of TTC35 recombinant protein (Accession Number: NM_014673). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116382(TTC35 Antibody) at dilution of 1:600



IP Result of anti-TTC35 (IP:Catalog No:116382, 3ug; Detection:Catalog No:116382 1:600) with mouse heart tissue lysate 4000ug.