

TTC7B antibody



Catalog Number: 116388

Product name

TTC7B antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TTC7B Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 50-53 kDa

Preparation

This antibody was obtained by immunization of TTC7B recombinant protein (Accession Number: NM_001010854). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:116388(TTC7B Antibody) at dilution of 1:600