

TROAP antibody



Catalog Number: 116402

Product name

TROAP antibody

Specificity

Human; other species not tested.

Antibody description

TROAP Rabbit Polyclonal antibody. Positive IF detected in MCF-7 cells. Positive WB detected in HeLa cells, PC-3 cells. Observed molecular weight by Western-blot: 110 kDa

Preparation

This antibody was obtained by immunization of TROAP recombinant protein (Accession Number: NM_005480). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:10-1:100

Validations



Immunofluorescent analysis of MCF-7 cells, using TROAP antibody Catalog No:116402 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116402(TROAP Antibody) at dilution of 1:300