

TRPM6 antibody

Catalog Number: 116412

Product name

TRPM6 antibody

Specificity

Human; other species not tested.

Antibody description

TRPM6 Rabbit Polyclonal antibody. Positive IF detected in HEK-293 cells.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_017662).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IF

Dilutions

Recommended Dilution:

IF: 1:10-1:100

Validations



Immunofluorescent analysis of HEK-293 cells using Catalog No:116412(TRPM6 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG