

# TSPAN13 antibody



Catalog Number: 116444

## Product name

TSPAN13 antibody

## Immunogen

[Human TSPAN13 Recombinant protein \(His tag\)](#)

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

TSPAN13 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue. Observed molecular weight by Western-blot: 28-35kd

## Preparation

This antibody was obtained by immunization of TSPAN13 recombinant protein (Accession Number: NM\_014399). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:116444(TSPAN13 antibody) at dilution of 1:800