

TSPAN13 antibody



Catalog Number: 116445

Product name

TSPAN13 antibody

Immunogen

[Human TSPAN13 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TSPAN13 Rabbit Polyclonal antibody. Positive IP detected in mouse lung tissue. Positive WB detected in mouse lung tissue. Observed molecular weight by Western-blot: 28-35kd

Preparation

This antibody was obtained by immunization of TSPAN13 recombinant protein (Accession Number: NM_014399). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:

WB: 1:200-1:2000


IP: 1:200-1:1000

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:116445(TSPAN13 antibody) at dilution of

1:300



IP Result of anti-TSPAN13 (IP:Catalog No:116445, 3ug; Detection:Catalog No:116445 1:300) with mouse lung tissue lysate 4000ug.