

# TSPAN13 antibody



Catalog Number: 116445

## Product name

TSPAN13 antibody

## Immunogen

[Human TSPAN13 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TSPAN13 Rabbit Polyclonal antibody. Positive IP detected in mouse lung tissue. Positive WB detected in mouse lung tissue. Observed molecular weight by Western-blot: 28-35kd

## Preparation

This antibody was obtained by immunization of TSPAN13 recombinant protein (Accession Number: NM\_014399). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:116445(TSPAN13 antibody) at dilution of



1:300

IP Result of anti-TSPAN13 (IP:Catalog No:116445, 3ug; Detection:Catalog No:116445 1:300) with mouse lung tissue lysate 4000ug.