

# TXNDC17 antibody



Catalog Number: 116449

## Product name

TXNDC17 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

TXNDC17 Rabbit Polyclonal antibody. Positive IF detected in Hela. Positive WB detected in Jurkat cells, HEK-293 cells, HeLa cells, human kidney tissue, K-562 cells, mouse pancreas tissue, human uterus tissue, L02 cells. Observed molecular weight by Western-blot: 14kd

## Preparation

This antibody was obtained by immunization of TXNDC17 recombinant protein (Accession Number: NM\_032731). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116449(TXNDC17 antibody) at dilution of



Immunofluorescent analysis of Hela cells, using TXNDC17 antibody Catalog No:116449 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).