

TXNDC3 antibody



Catalog Number: 116450

Product name

TXNDC3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TXNDC3 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, human testis tissue. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 67 kDa

Preparation

This antibody was obtained by immunization of TXNDC3 recombinant protein (Accession Number: NM_016616). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116450(TXNDC3 antibody) at dilution of



1:1000

IP Result of anti-TXNDC3 (IP:Catalog No:116450, 4ug; Detection:Catalog No:116450 1:1000) with mouse testis tissue lysate 4000ug.