

UACA antibody



Catalog Number: 116464

Product name

UACA antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

UACA Rabbit Polyclonal antibody. Positive IF detected in A375 cells. Positive WB detected in HeLa cells, A375 cells, HEK-293 cells, Neuro-2a cells, Y79 cells. Observed molecular weight by Western-blot: 163 kDa

Preparation

This antibody was obtained by immunization of UACA recombinant protein (Accession Number: NM_018003). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116464(UACA Antibody) at dilution of 1:2000



Immunofluorescent analysis of A375 cells using Catalog No:116464(UACA Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)