

UAP56 antibody



Catalog Number: 116466

Product name

UAP56 antibody

Recommended Dilution:

WB: 1:500-1:5000

Immunogen

IP: 1:200-1:2000

[Human BAT1 Recombinant protein \(GST tag\)](#)

IHC: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

UAP56 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive IP detected in Jurkat cells. Positive WB detected in Jurkat cells, rat spleen tissue. Observed molecular weight by Western-blot: 48 kDa




Preparation

This antibody was obtained by immunization of UAP56 recombinant protein (Accession Number: NM_004640). Purification method: Antigen affinity purified.

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116466(BAT1 antibody) at dilution of 1:1000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:116466(BAT1 antibody) at dilution of 1:100 (under 10x lens)

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IP


Dilutions

UAP56 antibody

Catalog Number: 116466



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:116466(BAT1 antibody) at dilution of 1:100



IP Result of anti-BAT1 (IP:Catalog No:116466, 4ug; Detection:Catalog No:116466 1:500) with Jurkat cells lysate 1200ug.