

# UAP56 antibody



Catalog Number: 116466

## Product name

UAP56 antibody

## Immunogen

[Human BAT1 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

UAP56 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive IP detected in Jurkat cells. Positive WB detected in Jurkat cells, rat spleen tissue. Observed molecular weight by Western-blot: 48 kDa

## Preparation

This antibody was obtained by immunization of UAP56 recombinant protein (Accession Number: NM\_004640). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


IP: 1:200-1:2000

IHC: 1:20-1:200

## Validations




Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116466(BAT1 antibody) at dilution of 1:1000




Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:116466(BAT1 antibody) at dilution of 1:100 (under 10x lens)

# UAP56 antibody

Catalog Number: 116466



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:116466(BAT1 antibody) at dilution of 1:100



(under 40x lens)

IP Result of anti-BAT1 (IP:Catalog No:116466, 4ug; Detection:Catalog No:116466 1:500) with Jurkat cells lysate 1200ug.