

UBASH3A-Specific antibody



Catalog Number: 116474

Product name

UBASH3A-Specific antibody

Recommended Dilution:

WB: 1:200-1:2000

Specificity

IP: 1:200-1:2000

IF: 1:10-1:100

Human, Mouse, Rat; other species not tested.

Antibody description

UBASH3A-Specific Rabbit Polyclonal antibody.
Positive WB detected in human spleen tissue.
Positive IP detected in mouse spleen tissue.
Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 69 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_018961).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

ELISA, WB, IF, IP


Dilutions

Validations



human spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:116474(UBASH3A-Specific antibody) at dilution of 1:400

Immunofluorescent analysis of HepG2 cells, using UBASH3A antibody Catalog No:116474 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

IP Result of anti-UBASH3A-Specific (IP:Catalog No:116474, 4ug; Detection:Catalog No:116474 1:500) with mouse spleen tissue lysate 4000ug.