

# UBASH3A-Specific antibody



Catalog Number: 116474

## Product name

UBASH3A-Specific antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

UBASH3A-Specific Rabbit Polyclonal antibody. Positive WB detected in human spleen tissue. Positive IP detected in mouse spleen tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 69 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_018961). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF, IP

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

IP: 1:200-1:2000


IF: 1:10-1:100

## Validations



human spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:116474(UBASH3A-Specific antibody) at dilution of 1:400

Immunofluorescent analysis of HepG2 cells, using UBASH3A antibody Catalog No:116474 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

IP Result of anti-UBASH3A-Specific (IP:Catalog No:116474, 4ug; Detection:Catalog No:116474 1:500) with mouse spleen tissue lysate 4000ug.