

UBASH3B/STS 1 antibody



Catalog Number: 116475

Product name

UBASH3B/STS 1 antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

UBASH3B/STS 1 Rabbit Polyclonal antibody. Positive IP detected in mouse kidney tissue. Positive WB detected in HEK-293 cells, mouse lung tissue. Positive IHC detected in human hepatocirrhosis tissue, human spleen tissue. Observed molecular weight by Western-blot: 66-70 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_032873). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC, IP, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog No:116475(UBASH3B antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog No:116475(UBASH3B antibody) at dilution of



1:100 (under 40x lens)

UBASH3B/STS 1 antibody

Catalog Number: 116475

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116475(UBASH3B Antibody) at dilution of

IP Result of anti-UBASH3B (IP:Catalog No:116475, 3ug; Detection:Catalog No:116475 1:300) with mouse kidney tissue lysate 6000ug.

