

# TTRAP antibody



Catalog Number: 116482

## Product name

TTRAP antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TTRAP Rabbit Polyclonal antibody. Positive WB detected in human spleen tissue. Positive IP detected in mouse spleen tissue. Observed molecular weight by Western-blot: 49kd

## Preparation

This antibody was obtained by immunization of TTRAP recombinant protein (Accession Number: NM\_016614). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



human spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:116482(TTRAP antibody) at dilution of 1:500



IP Result of anti-TTRAP (IP:Catalog No:116482, 4ug; Detection:Catalog No:116482 1:500) with mouse spleen tissue lysate 4800ug.