

# TULP2 antibody



Catalog Number: 116492

---

## Product name

TULP2 antibody

## Specificity

Human; other species not tested.

## Antibody description

TULP2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, human liver tissue, human lung tissue, Jurkat cells. Observed molecular weight by Western-blot: 59kd

## Preparation

This antibody was obtained by immunization of TULP2 recombinant protein (Accession Number: NM\_003323). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: N/A

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116492(TULP2 antibody) at dilution of 1:300