

# TULP3 antibody



Catalog Number: 116493

## Product name

TULP3 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200


## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TULP3 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells. Positive IP detected in SH-SY5Y cells. Positive IHC detected in human ovary tumor tissue. Observed molecular weight by Western-blot: 60 kDa

## Validations




## Preparation

This antibody was obtained by immunization of TULP3 recombinant protein (Accession Number: NM\_003324). Purification method: Antigen affinity purified.

SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:116493(TULP3 antibody) at dilution of 1:800

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT


## Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:116493(TULP3 antibody) at dilution of 1:100 (under 10x lens)

## Ig Type

Rabbit IgG



## Applications

ELISA, WB, IHC, IP

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Immunohistochemical of paraffin-embedded

# TULP3 antibody

Catalog Number: 116493

---

human ovary tumor using Catalog  
No:116493(TULP3 antibody) at dilution of 1:100

IP Result of anti-TULP3 (IP:Catalog No:116493,  
4ug; Detection:Catalog No:116493 1:800) with  
SH-SY5Y cells lysate 1600ug.

