

# TULP4 antibody



Catalog Number: 116494

## Product name

TULP4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TULP4 Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue, human brain tissue. Positive IHC detected in human kidney tissue. Observed molecular weight by Western-blot: 169kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_020245). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:116494(TULP4 antibody) at dilution of 1:800