

**Product name**

Twinkle antibody

**Specificity**

Human; other species not tested.

**Antibody description**

Twinkle Rabbit Polyclonal antibody. Positive IHC detected in human pancreas tissue, human skeletal muscle tissue. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 60kd

**Preparation**

This antibody was obtained by immunization of Twinkle recombinant protein (Accession Number: BC033762). Purification method: Antigen affinity purified.

**Formulation**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

**Storage**

Store at -20°C. DO NOT ALIQUOT

**Clonality**

Polyclonal

**Ig Type**

Rabbit IgG

**Applications**


ELISA, WB, IHC

**Dilutions**


Recommended Dilution:

WB: 1:200-1:2000

IHC: 1:20-1:200

**Validations**

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116500(PEO1 antibody) at dilution of 1:800



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using Catalog No:116500(PEO1 Antibody) at dilution of 1:50