

TWIST1-specific antibody



Catalog Number: 116502

Product name

TWIST1-specific antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TWIST1-specific Rabbit Polyclonal antibody.
Positive IP detected in MDA-MB-453s cells.
Positive WB detected in MDA-MB-453s cells, A549 cells, mouse brain tissue, Raji cells. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of TWIST1-specific recombinant protein (Accession Number: XM_011515496). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



MDA-MB-453s cells were subjected to SDS PAGE followed by western blot with Catalog No:116502(TWIST1 Antibody) at dilution of 1:600



IP Result of anti-TWIST1 (IP:Catalog No:116502, 4ug; Detection:Catalog No:116502 1:500) with MDA-MB-453s cells lysate 920ug.