

TWIST2 antibody



Catalog Number: 116503

Product name

TWIST2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TWIST2 Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells, MCF7 cells, mouse liver tissue. Observed molecular weight by Western-blot: 21kDa, 60 kDa

Preparation

This antibody was obtained by immunization of TWIST2 recombinant protein (Accession Number: NM_001271893). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:116503(TWIST2 antibody) at dilution of 1:500