

TXNDC11 antibody



Catalog Number: 116506

Product name

TXNDC11 antibody

Specificity

Human; other species not tested.

Antibody description

TXNDC11 Rabbit Polyclonal antibody. Positive WB detected in DU 145 cells. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of TXNDC11 recombinant protein (Accession Number: BC013727). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



DU 145 cells were subjected to SDS PAGE followed by western blot with Catalog No:116506(TXNDC11 Antibody) at dilution of 1:1000