

UBE1L antibody



Catalog Number: 116514

Product name

UBE1L antibody

IHC: 1:20-1:200

Specificity

Human; other species not tested.

Antibody description

UBE1L Rabbit Polyclonal antibody. Positive WB detected in human liver tissue. Positive IHC detected in human endometrial cancer tissue. Observed molecular weight by Western-blot: 112kd

Preparation

This antibody was obtained by immunization of UBE1L recombinant protein (Accession Number: NM_003335). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:116514(UBA7 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human endometrial cancer using Catalog No:116514(UBA7 antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemical of paraffin-embedded human endometrial cancer using Catalog No:116514(UBA7 antibody) at dilution of 1:50

UBE1L antibody



Catalog Number: 116514

(under 10x lens)