

# UBE2A antibody



Catalog Number: 116515

## Product name

UBE2A antibody

ELISA, WB, IHC

## Specificity

Human, Mouse, Rat; other species not tested.

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IHC: 1:20-1:200

## Antibody description

UBE2A Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue. Positive WB detected in HEK-293 cells, COLO 320 cells, human kidney tissue, human placenta tissue. Observed molecular weight by Western-blot: 17 kDa

## Preparation

This antibody was obtained by immunization of UBE2A recombinant protein (Accession Number: NM\_003336). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality


Polyclonal

## Ig Type


Rabbit IgG

## Applications

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116515(UBE2A antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:116515(UBE2A antibody) at dilution of 1:50 (under 10x lens)