

UBE2E2 antibody

Catalog Number: 116522

Product name

UBE2E2 antibody

IHC: 1:20-1:200

IF: 1:10-1:100


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description


UBE2E2 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue. Positive IF detected in HepG2 cells. Positive WB detected in human thyroid gland tissue, human brain tissue, mouse embryo tissue. Observed molecular weight by Western-blot: 22 kDa



Preparation

This antibody was obtained by immunization of UBE2E2 recombinant protein (Accession Number: NM_001370225). Purification method: Antigen affinity purified.

human thyroid gland tissue were subjected to SDS PAGE followed by western blot with Catalog No:116522(UBE2E2 antibody) at dilution of 1:300



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Immunohistochemical of paraffin-embedded human gliomas using Catalog No:116522(UBE2E2 antibody) at dilution of 1:50 (under 10x lens)

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:

WB: 1:200-1:2000



Immunofluorescent analysis of HepG2 cells, using UBE2E2 antibody Catalog No:116522 at 1:25

UBE2E2 antibody



Catalog Number: 116522

dilution and Rhodamine-labeled goat anti-rabbit IgG (red).