

UCP5 antibody



Catalog Number: 116548

Product name

Polyclonal

UCP5 antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

UCP5 Rabbit Polyclonal antibody. Positive WB detected in fetal human brain tissue, PC-3 cells. Observed molecular weight by Western-blot: 36-40 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Preparation

This antibody was obtained by immunization of UCP5 recombinant protein (Accession Number: NM_003951). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116548(SLC25A14 Antibody) at dilution of 1:300