

**Product name**

UCRC antibody

**Specificity**

Human, Mouse, Rat; other species not tested.

**Antibody description**

UCRC Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, human liver tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 7 kDa

**Preparation**

This antibody was obtained by immunization of UCRC recombinant protein (Accession Number: NM\_013387). Purification method: Antigen affinity purified.

**Formulation**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

**Storage**

Store at -20°C. DO NOT ALIQUOT

**Clonality**

Polyclonal

**Ig Type**

Rabbit IgG


**Applications**

ELISA, WB

**Dilutions**

Recommended Dilution:

WB: 1:500-1:5000

**Validations**

human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:116549(UCRC antibody) at dilution of 1:300