

# UGCGL2 antibody



Catalog Number: 116555

## Product name

Polyclonal

UGCGL2 antibody

## Ig Type

## Specificity

Rabbit IgG

human; other species not tested.

## Applications

## Antibody description

ELISA, WB

UGCGL2 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 170kd

## Dilutions

Recommended Dilution:

## Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of UGCGL2 recombinant protein (Accession Number: BC032302). Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116555(UGCGL2 antibody) at dilution of 1:300