

UHMK1 antibody

Catalog Number: 116562

Product name

UHMK1 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

UHMK1 Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer tissue. Positive WB detected in mouse thymus tissue, mouse brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 47kd

Preparation

This antibody was obtained by immunization of UHMK1 recombinant protein (Accession Number: NM_144624). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:116562(UHMK1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:116562(UHMK1 antibody) at dilution of 1:50 (under 10x lens)



IP Result of anti-UHMK1 (IP:Catalog No:116562,

UHMK1 antibody



Catalog Number: 116562

4ug; Detection:Catalog No:116562 1:500) with mouse brain tissue lysate 2640ug.