

# UMP/CMP kinase antibody



Catalog Number: 116569

## Product name

UMP/CMP kinase antibody

ELISA, WB, IHC

## Specificity

Human, Mouse, Rat; other species not tested.

## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Antibody description

UMP/CMP kinase Rabbit Polyclonal antibody.  
Positive WB detected in human pancreas tissue,  
human skeletal muscle tissue. Positive IHC  
detected in human pancreas cancer tissue.  
Observed molecular weight by Western-blot: 26  
kDa

## Validations




## Preparation

This antibody was obtained by immunization of  
UMP/CMP kinase recombinant protein (Accession  
Number: NM\_016308). Purification method:  
Antigen affinity purified.

human pancreas tissue were subjected to SDS  
PAGE followed by western blot with Catalog  
No:116569(CMPK1 antibody) at dilution of 1:300

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH  
7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

Immunohistochemical of paraffin-embedded  
human pancreas cancer using Catalog  
No:116569(CMPK1 antibody) at dilution of 1:100  
(under 10x lens)

## Applications