

USH1G antibody

Catalog Number: 116582

Product name

USH1G antibody

Specificity

Human; other species not tested.

Antibody description

USH1G Rabbit Polyclonal antibody. Positive IP detected in NIH/3T3 cells. Observed molecular weight by Western-blot: 52-55 kDa

Preparation

This antibody was obtained by immunization of USH1G recombinant protein (Accession Number: NM_173477). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IP

Dilutions

Recommended Dilution:

IP: 1:200-1:2000

Validations



IP Result of anti-USH1G (IP:Catalog No:116582, 4ug; Detection:Catalog No:116582 1:500) with NIH/3T3 cells lysate 1200ug.