

USP15 antibody

Catalog Number: 116592

Product name

USP15 antibody

Immunogen

[Human USP15 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

USP15 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue, HEK-293 cells, HeLa cells. Positive IHC detected in human breast cancer tissue, human ovary tumor tissue. Observed molecular weight by Western-blot: 112 kDa

Preparation

This antibody was obtained by immunization of USP15 recombinant protein (Accession Number: NM_001252079). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP, IHC

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


IP: 1:500-1:5000

IHC: 1:50-1:500

Validations




mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116592(USP15 Antibody) at dilution of 1:1000




IP Result of anti-USP15 (IP:Catalog No:116592, 4ug; Detection:Catalog No:116592 1:1000) with mouse brain tissue lysate 4000ug.

USP15 antibody

Catalog Number: 116592



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:116592(USP15 Antibody) at dilution of 1:200 (under 10x lens). 

Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:116592(USP15 Antibody) at dilution of 1:200 (under 40x lens).