

USP18 antibody

Catalog Number: 116594

Product name

USP18 antibody

Specificity

Human; other species not tested.

Antibody description

USP18 Rabbit Polyclonal antibody. Positive WB detected in MCF-7. Observed molecular weight by Western-blot: 43 kDa, 70 kDa

Preparation

This antibody was obtained by immunization of USP18 recombinant protein (Accession Number: NM_017414). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: N/A

Validations



WB result of USP18 antibody(Catalog No:116594) with MCF-7 cells (Control and CRISPR USP18).